

# SILPAK, Inc

470 E. BONITA AVE., POMONA, CA 91767  
PH (909) 625-0056 [WWW.SILPAK.COM](http://WWW.SILPAK.COM) FX (909) 625-0082

## R-2000 V (Platinum RTV Retarder) Product Data Sheet

R 2000 V is additive used in Platinum RTVs (R-2000 Series) to extend and control material Gel Time. Add 1% to Part B (Catalyst). Amount to add is calculated on Total Weight of Part B.

Available Sizes: 2oz, 8oz, 1 lb Bottles

### MEASUREMENTS

- 1% = 4.5 gram/lb
- 2% = 9 grams/lb

**Note:** Gel Time (Pot Life) will vary with each of the *R-2000 Series Platinum RTV* systems. Use the above recommended levels as a general starting point. Please test workability of silicone to assure good results.  
*\*\*Environmental conditions such as ambient temperature, especially hot weather, will also shorten the amount of Gel and Cure Times.*

THE INFORMATION AND DATA CONTAINED HEREIN ARE BASED ON INFORMATION WE BELIEVE RELIABLE. EACH USER OF THE MATERIAL SHOULD THOROUGHLY TEST ANY APPLICATION AND INDEPENDENTLY CONCLUDE SATISFACTORY PERFORMANCE BEFORE COMMERCIALIZING. SUGGESTIONS OF USES SHOULD NOT BE TAKEN AS INDUCEMENTS TO INFRINGE ON ANY PARTICULAR PATENT.